


PCR Screening for TLR4 Targeted Mice


Purpose: To identify Toll receptor 4 (TLR4) deficient mice.

Principle: The murine Tlr4 gene contains 3 exons (transcript length is 3843 bps, and translation length is 835 bps) and is located in mouse chromosome 4 (ensemble: ENSMUSG00000039005, NCBI: NM_021297.2, location: 66,313,972-66,328,954). The gene has been disrupted by insertion of a Neomycin (**Neo**) cassette (unknown origin) sequence in replacement of exon 3. The screen amplifies a section of intron 2-3 (Forward primer HC-26 attaches in intron 2-3 from position 4372-4399) and either the **Neo** cassette to yield a band specific to the KO (1300bp) or a section of exon 3 (Reverse primer HC-27 attaches in exon 3 from position 503-477) for the wild-type band (1300bp).

Wild-Type Band (~1300bp)



Targeted Band (~1300bp)




Primers:

TLR4 HC-26 F: 5' -TGTTGCCCTTCAGTCACAGAGACTCTG-3' 27-mer

TLR4 HC-27 R: 5' -CGTGTAACCAGCCAGGTTTTGAAGGC-3' 27-mer

TLR4 HC-28 NEO: 5' -ATCGCCTTCTATCGCCTTCTTGACGAG-3' 27-mer



Wild type gel is run on the top,
deficient gel is run on the bottom.
All mice in this screen were TLR4+/+.

Wild-type locus will produce a 1300 bp product
Targeted locus will produce a 1300 bp product

Reaction Components:

Wild-type Reaction

DNA template	2.0 μ l
TLR4 HC-26 Forward	0.2 μ l
TLR4 HC-27 Reverse	0.2 μ l
*Promega Master Mix	10.0 μ l
*Taq	0.1 μ l
H ₂ O	<u>7.5 μl</u>
	20.0 μ l total

PCR Program:

94° for 4 minutes	
94° for 1 minute	} 30 cycles
62° for 1 minute	
72° for 2 minutes	
72° for 10 Minutes	

Targeted Reaction

DNA template	2.0 μ l
TLR4 HC-26 Forward	0.2 μ l
TLR4 HC-28 NEO	0.2 μ l
*Promega Master Mix	10.0 μ l
*Taq	0.1 μ l
H ₂ O	<u>7.5 μl</u>
	20.0 μ l total

94° for 4 minutes	
94° for 1 minute	} 30 cycles
62° for 1 minute	
72° for 2 minutes	
72° for 10 Minutes	

*Note: If you are using Promega GoTaq Green Master Mix, you do not have to add Taq polymerase

Electrophoresis Conditions:

Since the two products come out at exactly the same (1300bp), you must make sure to run on a two row gel, or make two gels. Run with 2% agarose loading 15-20ul at 100-110 Volts for ~30 minutes in order to achieve separation.

This protocol was adapted from Dr. Peter Tobias at the Scripps Institute.

Modified: 09-16-07 APO

Updated: 07-10-08 APO